



Plasmodium falciparum PF10Nter Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1139
ORIGIN SPECIES	Plasmodium falciparum
FRAGMENT	Partial
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	95,44 kDa
FORM	liquid
BUFFER	TrisHC 50mMl, 10mM reduced glutathion pH8
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q8IK14
ALIAS/SYNONYMS	PF10Nter , Tryptophan-rich antigen 3, putative, PF10_0026, PF3D7_1002200

