



## Plasmodium falciparum PFL1260 Recombinant Protein

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-P1140
<b>ORIGIN SPECIES</b>	Plasmodium falciparum
<b>FRAGMENT</b>	Full-length
<b>TAG</b>	Yes
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Prokaryotic expression
<b>HOST</b>	Escherichia coli (E.coli)
<b>MOLECULAR WEIGHT</b>	37,55 kDa
<b>FORM</b>	liquid
<b>BUFFER</b>	100mM Na <sub>2</sub> HPO <sub>4</sub> , TrisHCl 10mM, Urea 8M
<b>STABILITY &amp; STORAGE</b>	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
<b>UNIPROT</b>	Q8I5F6
<b>ALIAS/SYNONYMS</b>	PFL1260, hydrolase / phosphatase, putative

