



Trypanosoma PLB Recombinant Protein N-Ter His tag

OVERVIEW

PRODUCT REFERENCE	PX-P1141
ORIGIN SPECIES	Trypanosoma
FRAGMENT	Full-length
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	33.31 kDa
FORM	liquid
BUFFER	PBS, Urea 8M
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	F9WF01
ALIAS/SYNONYMS	PLB, unnamed protein product

