



Drosophila RbF1 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1145
ORIGIN SPECIES	Drosophila
FRAGMENT	Partial
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	94,87 kDa
FORM	liquid
BUFFER	TrisHC 50mMl, 10mM glutathione reduced pH8
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q24472
ALIAS/SYNONYMS	RbF1, Retinoblastoma-family protein, CG7413, rb, pRb, rbf, dRBF, RB, Rb1, Dmel\CG7413, EG:34F3.3, Rb, RBF, RbF, rbf1, CG7413, FBF

