



Ciona intestinalis SLC Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1150
ORIGIN SPECIES	Ciona intestinalis
FRAGMENT	Partial
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	12,97 kDa
FORM	liquid
BUFFER	PBS, imidazole 200mM if native conditions. PBS, imidazole 10mM, Urea 5M if denaturing conditions
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	T2DR23
ALIAS/SYNONYMS	SLC, sulfate transporter Ci-Slc26a alpha

