



Tomato SI_eIF4E2 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1151
ORIGIN SPECIES	Tomato
FRAGMENT	Partial
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	23,10 kDa
FORM	liquid
BUFFER	PBS, Urea 8M
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	XP_004231545.1
ALIAS/SYNONYMS	SI_eIF4E2, eukaryotic translation initiation factor 4E-1-like

