



S. pombe Smc3 (Nter) Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1153
ORIGIN SPECIES	S. pombe
FRAGMENT	Partial
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	74,38 kDa
FORM	liquid
BUFFER	PBS, 6M Urea
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	O42649
ALIAS/SYNONYMS	Smc3 (Nter), mitotic cohesin complex subunit Psm3

