

Human SR-BI Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P1156
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Partial
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	23,24 kDa
FORM	liquid
BUFFER	PBS, Urea 8M, imidazole 400mM
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q8WTV0
ALIAS/SYNONYMS	SR-BI [310-882], scavenger receptor class B, member 1

