



## Rhesus monkey TAFA2 Recombinant Protein

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-P1158
<b>ORIGIN SPECIES</b>	Rhesus monkey
<b>FRAGMENT</b>	Partial
<b>TAG</b>	Yes
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Prokaryotic expression
<b>HOST</b>	Escherichia coli (E.coli)
<b>MOLECULAR WEIGHT</b>	13,16 kDa
<b>FORM</b>	liquid
<b>BUFFER</b>	PBS, imidazole 300mM in natives conditions. PBS, imidazole 10mM, Urea 8M in denaturing conditions
<b>STABILITY &amp; STORAGE</b>	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
<b>UNIPROT</b>	G7N7I3
<b>ALIAS/SYNONYMS</b>	TAFA2, Chemokine-like protein TAFA-2

