



Human UHRF1(400-660) tagged Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P2053
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Partial
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	30,19 kDa
FORM	liquid
BUFFER	100mM Na ₂ HPO ₄ , TrisHCl 10mM, Urea 8M in denaturing conditions. PBS if renaturation.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q96T88
ALIAS/SYNONYMS	ICBP[400-660], transcription factor ICBP90

