



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Human Chromogranin A Partial [19-457] Recombinant Protein C- His Tag

OVERVIEW

PRODUCT REFERENCE	PX-P2060
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Partial
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	68.16KDa
FORM	Lyophilized
BUFFER	Tris 50mM and 500mM NaCl pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P10645
ALIAS/SYNONYMS	ChromograninA Partial [19-457], chromogranin A (parathyroid secretory protein 1), CHGA, Pituitary secretory protein I, SP-I,



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.