



Staphylococcus AmiR1R2 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P2063
ORIGIN SPECIES	Staphylococcus
FRAGMENT	Unknown
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	40.59 kDa
FORM	Frozen
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	KEK49756.1
ALIAS/SYNONYMS	AmiR1R2, Bifunctional N-acetylmuramoyl-L-alanine amidase/endo-beta-N-acetylglucosaminidase, Atl

