



Drosophila CG30417 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P2066
ORIGIN SPECIES	Drosophila
FRAGMENT	Partial
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	31.95kDa
FORM	Frozen
BUFFER	PBSurea +8M, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9W1S2
ALIAS/SYNONYMS	CG30417, CG30417 ou IP01291p

