



## Human eIF4G Recombinant Protein

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-P2069
<b>ORIGIN SPECIES</b>	Homo sapiens (Human)
<b>FRAGMENT</b>	Glu557~Lys646
<b>TAG</b>	N-ter His&Trx Tag
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Prokaryotic expression
<b>HOST</b>	Escherichia coli (E.coli)
<b>MOLECULAR WEIGHT</b>	23.49kDa
<b>FORM</b>	liquid
<b>BUFFER</b>	PBS, imidazole 200mM + 50% glycérol
<b>STABILITY &amp; STORAGE</b>	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
<b>UNIPROT</b>	EAW78268.1
<b>ALIAS/SYNONYMS</b>	EIF-4G1, EIF4F, EIF4G, EIF4GI, P220, PARK18

