

Caenorhabditis F52B5 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P2070
ORIGIN SPECIES	Caenorhabditis
FRAGMENT	Partial
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	35.69 kDa
FORM	Frozen
BUFFER	PBS, pH 7.5 urea +8M
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q20643
ALIAS/SYNONYMS	F52B5, F52B5.2

