



Human DCLK3 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P2074
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Partial
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	74.95kDa
FORM	Frozen
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9C098
ALIAS/SYNONYMS	Human DCLK3, doublecortin-like kinase 3, Serine/threonine-protein kinase DCLK3, Doublecortin domain-containing protein 3C. Doublecortin-like and CAM kinase-like 3, Doublecortin-like kinase 3

