



Plasmodium MSP1-p19 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P2080
ORIGIN SPECIES	Plasmodium
FRAGMENT	Partial
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	23.79 kDa
FORM	Lyophilized
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	X2J7U7
ALIAS/SYNONYMS	MSP1-p19, Merozoite surface protein 1

