



## Synthetic construction GST 26-1 Recombinant Protein

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-P2092
<b>ORIGIN SPECIES</b>	synthetic construction
<b>FRAGMENT</b>	
<b>TAG</b>	
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Prokaryotic expression
<b>HOST</b>	Escherichia coli (E.coli)
<b>MOLECULAR WEIGHT</b>	39.38kDa
<b>FORM</b>	liquid
<b>BUFFER</b>	
<b>STABILITY &amp; STORAGE</b>	
<b>UNIPROT</b>	P08515
<b>ALIAS/SYNONYMS</b>	GST 26-1 (Sta1)

