



Toxoplasma gondii TgRA15 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P2097
ORIGIN SPECIES	Toxoplasma gondii
FRAGMENT	
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	29.39 kDa
FORM	liquid
BUFFER	
STABILITY & STORAGE	
UNIPROT	S8EY82
ALIAS/SYNONYMS	TgRA15, RAP domain-containing protein

