



VSVG Protein- Drosophila VSVG NJ Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P2099
ORIGIN SPECIES	Drosophila
FRAGMENT	Full-length
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Insect
MOLECULAR WEIGHT	59.11kDa
FORM	Frozen
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P15425
ALIAS/SYNONYMS	VSVG, neither inactivation nor afterpotential A, Peptidyl-prolyl cis-trans isomerase, rhodopsin-specific isozyme, PPIase, Rotamase

