

Staphylococcus 6His-Protein A

OVERVIEW

CATALOG NO	PX-P2101
ORIGIN SPECIES	Staphylococcus
FRAGMENT	Partial
MOLECULAR WEIGHT	33.83 KDa
DELIVERED WITH TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
PURITY	90%
UNIPROT ID	P38507
FORM	Liquid
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)

More info

SYNONYMS

6His-Protein A

UNIPROT LINK

<http://www.uniprot.org/uni-prot/P38507>

SEQUENCE

MGAQHDEAQQNAFYQVLNMPNLNADQRNGFIQSLKDDPSQSANVLGEAQKLNSQAPKADAQQNKFNKDQQSAFYEILNMPNLNEEQRNGFIQSLKDDPSQSTNVLGEAKKLNESQAPKADNMFNKEQQNAFYEILNMPNLNEEQRNGFIQSLKDDPSQSANLLAEAKKLNESQAPKADNMFNKEQQNAFYEILHLPNLNEEQRNGFIQSLKDDPSQSANLLAEAKKLNDAPKADNMFNKEQQNAFYEILHLP