



Proteogenix

Antibodies and Proteins

Proteogenix SAS
Espace Européen de l'Entreprise
19 rue de La Haye
67300 Schiltigheim, France
Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88
contact@proteogenix.fr
www.Proteogenix.science

Streptococcus 6His-Protein G

OVERVIEW

CATALOG NO	PX-P2102
ORIGIN SPECIES	Streptococcus
FRAGMENT	Partial
MOLECULAR WEIGHT	20.96 KDa
DELIVERED WITH TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Procaryotic expression
HOST	Escherichia coli (E.coli)
PURITY	90%
UNIPROT ID	
FORM	Liquid
BUFFER	PBS,pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)

More info

SYNONYMS

6His-Protein G

UNIPROT LINK

SEQUENCE

MGKGETTTEAVDAATAEKVFKQYANDNGVDGEWTYDDATKTFTVTEKPEVIDASELTPAVTTYKLVINGKTLKGETTTEAVDAATAEKV
FKQYANDNGVDGEWTYDDATKTFTVTEKPEVIDASELTPAVTTYKLVINGKTLKGETTTKAVDAETAEKAFKQYANDNGVDGVWYDD
ATKTFTVTEHHHHHH