



Human Pepsin A-3(PGI) Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P3000
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Full-length
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	44.68 kDa
FORM	Frozen
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P0DJD8
ALIAS/SYNONYMS	EC 3.4.23; EC 3.4.23.1; FLJ58952; FLJ77962; pepsin A; pepsinogen 4, group I (pepsinogen A); Pepsinogen A; pepsinogen A4; PGA3; PGA5

