



## Human Sex hormone-binding globulin (SHBG) Recombinant Protein

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-P3004
<b>ORIGIN SPECIES</b>	Homo sapiens (Human)
<b>FRAGMENT</b>	Full-length
<b>TAG</b>	Yes
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Eukaryotic expression
<b>HOST</b>	Mammalian cells
<b>MOLECULAR WEIGHT</b>	44.82 kDa
<b>FORM</b>	Frozen
<b>BUFFER</b>	PBS, pH 7.5
<b>STABILITY &amp; STORAGE</b>	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
<b>UNIPROT</b>	P04278
<b>ALIAS/SYNONYMS</b>	ABP, TEBG, SBP, Sex steroid-binding protein, Testis-specific androgen-binding protein, Testosterone-estradiol-binding globulin, Testosterone-estrogen-binding globulin, Sex Hormone Binding Globulin, SHBG

