

Human IL6 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P3013
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Full-length
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	49.8 kDa
FORM	Frozen
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P05231
ALIAS/SYNONYMS	Interleukin 6, IL-6, MGI2-A, MGI2A, HGF, BSF2, HSF, IFNB2, B-Cell Stimulatory Factor-2, B-cell stimulatory factor 2, Hybridoma/Plasmacytoma Growth Factor proteins, Hepatocyte Stimulating Factor, Cytotoxic T-Cell Differentiation Factor, CTL differentiation

