

Dog Programmed death ligand 1(PDL1)-CD274 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P3014
ORIGIN SPECIES	Dog
FRAGMENT	Partial
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	25.94 kDa
FORM	Frozen
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	E2RKZ5
ALIAS/SYNONYMS	CD274, PDL1, B7-H, B7H1, PD-L1, PDCD1L1, PDCD1-LG1, Programed Death Ligand 1, Programmed Cell Death Protein 1 Ligand 1, PDCD1LG1