



FSH protein - Chinese sturgeon Follicle-stimulating hormone(Follicle) Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P3016
ORIGIN SPECIES	Chinese sturgeon
FRAGMENT	Full-length
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	39.85 kDa
FORM	Frozen
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	B2CKR3
ALIAS/SYNONYMS	Follicle Stimulating Hormone, FSH

