



## Human FGF1 Recombinant Protein

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-P3028
<b>ORIGIN SPECIES</b>	Homo sapiens (Human)
<b>FRAGMENT</b>	Full-length
<b>TAG</b>	Yes
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Prokaryotic expression
<b>HOST</b>	Escherichia coli (E.coli)
<b>MOLECULAR WEIGHT</b>	18.78 kDa
<b>FORM</b>	Frozen
<b>BUFFER</b>	PBS,pH 7.5
<b>STABILITY &amp; STORAGE</b>	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
<b>UNIPROT</b>	P05230
<b>ALIAS/SYNONYMS</b>	ECGF, AFGF, ECGFA, ECGFB, FGF-Alpha, FGFA, HBGF1, ECGFB, HBGF1, Acidic Fibroblast Growth Factor proteins, heparin-binding Growth Factor proteins 1, Endothelial Cell Growth Factor proteins, Beta

