



Human PDGFB Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P3029
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Partial
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	26.26 kDa
FORM	Frozen
BUFFER	PBS,pH 7.5 urea +8M
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P01127
ALIAS/SYNONYMS	PDGF-B, PDGF2, SIS, SSV, C-sis, Simian Sarcoma Viral(v-sis)Oncogene Homolog, Platelet Derived Growth Factor proteins Beta Polypeptide, Proto-oncogene c-Sis, Becaplermin, Platelet Derived Growth Factor proteins Subunit B,PDGFB

