



Human TSG101 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P3051
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Partial
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	17.44 kDa
FORM	Frozen
BUFFER	PBS,pH 7.5 urea +8M
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q99816
ALIAS/SYNONYMS	TSG10, VPS23, Tumor Susceptibility Gene 10, ESCRT-I complex subunit TSG101, Tumor Susceptibility Gene 101

