

## Human Eosinophil cationic protein(ECP) Recombinant Protein

### OVERVIEW

---

|                           |                      |
|---------------------------|----------------------|
| <b>PRODUCT REFERENCE</b>  | PX-P3057             |
| <b>ORIGIN SPECIES</b>     | Homo sapiens (Human) |
| <b>FRAGMENT</b>           | Full-length          |
| <b>TAG</b>                |                      |
| <b>DELIVERY CONDITION</b> | Dry Ice              |

### PRODUCT INFORMATION

---

|                                |                                                                                                                                                |
|--------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>EXPRESSION SYSTEM</b>       | Eukaryotic expression                                                                                                                          |
| <b>HOST</b>                    | Mammalian cells                                                                                                                                |
| <b>MOLECULAR WEIGHT</b>        | 44.53 kDa                                                                                                                                      |
| <b>FORM</b>                    | Frozen                                                                                                                                         |
| <b>BUFFER</b>                  | PBS,pH 7.5                                                                                                                                     |
| <b>STABILITY &amp; STORAGE</b> | 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection) |
| <b>UNIPROT</b>                 | P12724                                                                                                                                         |
| <b>ALIAS/SYNONYMS</b>          | Rnase-A3, RnaseA3, ECP, RNS3, Ribonuclease,RNase A Family 3, Eosinophil Cationic Protein                                                       |