

Human VEGFR2 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P3059
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Partial
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	86.24 kDa
FORM	Frozen
BUFFER	PBS,pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P35968
ALIAS/SYNONYMS	Vascular endothelial Growth Factor proteins receptor 2, KDR,VEGFR-2, Fetal liver kinase 1, FLK-1, Kinase insert domain receptor, Protein-tyrosine kinase receptor flk-1, CD309

