



Elizabethkingia miricola PNGaseF Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P3070
ORIGIN SPECIES	Elizabethkingia miricola
FRAGMENT	Full-length
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	34.78 kDa
FORM	Frozen
BUFFER	Tris-HCl 20mM, NaCl 50mM, EDTA 5mM, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P21163
ALIAS/SYNONYMS	Peptide -N-Glycosidase F, PNGase F, PNGaseF

