



## Human T-cell immunoreceptor with Ig and ITIM domains (TIGIT) recombinant protein

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-P4036
<b>ORIGIN SPECIES</b>	Homo sapiens (Human)
<b>FRAGMENT</b>	Met1-Pro141
<b>TAG</b>	C-Terminal His Tag
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Eukaryotic expression
<b>HOST</b>	Mammalian cells
<b>MOLECULAR WEIGHT</b>	26.32 kDa
<b>FORM</b>	Lyophilized
<b>BUFFER</b>	PBS pH 7.4
<b>STABILITY &amp; STORAGE</b>	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
<b>UNIPROT</b>	Q495A1
<b>ALIAS/SYNONYMS</b>	TIGIT, VSIG9, VSTM3, V-Set And Immunoglobulin Domain Containing 9, V-Set And Transmembrane Domain Containing 3

