



## Human PDL1, B7-H1, CD274 recombinant protein

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-P4038
<b>ORIGIN SPECIES</b>	Homo sapiens (Human)
<b>FRAGMENT</b>	Partial
<b>TAG</b>	C-Terminal His Tag
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Eukaryotic expression
<b>HOST</b>	Mammalian cells
<b>MOLECULAR WEIGHT</b>	28.73kDa
<b>FORM</b>	Lyophilized
<b>BUFFER</b>	50 mM Tris-HCl pH 8, 150 mM NaCl
<b>STABILITY &amp; STORAGE</b>	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
<b>UNIPROT</b>	Q9NZQ7
<b>ALIAS/SYNONYMS</b>	PDL1, Programmed cell death 1 ligand 1, CD274, B7-H, B7H1, PD-L1, PDCD1L1, PDCD1-LG1, Programed Death Ligand 1

