



Human Calcitonin-Procalcitonin (PCT) recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P4042
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Full-length
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	16.4kDa
FORM	Lyophilized
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P01258
ALIAS/SYNONYMS	PCT, CALCA, Alpha CGRP, alpha type CGRP, CALC 1, CALC A, CALC1, calcitonin 1, calcitonin gene related peptide 1, calcitonin gene related peptide I, calcitonin related polypeptide alpha, CGRP, CGRP 1, CGRP I, CGRP1, CT, katalcalcin, KC

