



Human CD28 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P4044
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Partial
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	42.63kDa
FORM	Lyophilized
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P10747
ALIAS/SYNONYMS	CD28, T-cell-specific surface glycoprotein CD28, Tp44, T-Cell-Specific Surface Glycoprotein, MGC138290

