



Proteogenix

Antibodies and Proteins

Proteogenix SAS
Espace Européen de l'Entreprise
19 rue de La Haye
67300 Schiltigheim, France
Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88
contact@proteogenix.fr
www.Proteogenix.science

Human CD28 recombinant protein

OVERVIEW

CATALOG NO	PX-P4044
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Partial
MOLECULAR WEIGHT	42.63kDa
DELIVERED WITH TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
PURITY	70%
UNIPROT ID	P10747
FORM	Liquid
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)

More info

SYNONYMS

CD28, T-cell-specific surface glycoprotein CD28, Tp44, T-Cell-Specific Surface Glycoprotein, MGC138290

UNIPROT LINK

<http://www.uniprot.org/uni/prot/P10747>

SEQUENCE

MLRLLLALNLFPSIQVTGNKILVKQSPMLVAYDNAVNLSCKYSYNLFSREFRASLHKGLDSAVEVCVVYGNYSQQLQVYSKTGFNC DGKL
GNESVTFYLQNLVYVNTDIYFCKIEVMYPPPYLDNEKSNGTIIHVKGKHLCPSPFLPGPSKPKDKTHTCPPCPAPPELLGGPSVFLFPPKPKDQL
MISRTPEVTCVVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIEK TISKAK
GQPREPQVYITLPPSREEMTKNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKTTTPVLDSDGSFFLYSKLTVDKSRWQQGNV FSCSVLHEA
LHNHYTQKSLSLSPGK