



Recombinant human VEGF Receptor 1 protein-FLT 1

OVERVIEW

PRODUCT REFERENCE	PX-P4058
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Partial
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	37.7kDa
FORM	Lyophilized
BUFFER	PBS, pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P17948
ALIAS/SYNONYMS	EC 2.7.10.1, FLT, FLT 1, Flt-1, FLT1, Fms like tyrosine kinase 1, Fms related tyrosine kinase 1, Fms related tyrosine kinase 1 (vascular endothelial Growth Factor proteins/vascular permeability factor receptor), Fms related tyrosine kinase 1 vascular endo

