



Mouse PAD4 protein-arginine deiminase type-4 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P4068
ORIGIN SPECIES	House Mouse
FRAGMENT	Full-length
TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Insect
MOLECULAR WEIGHT	77.52kDa
FORM	Frozen
BUFFER	50mM Tris-Cl pH 8.0, 150mM NaCl
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9Z183
ALIAS/SYNONYMS	Pad4, Pdi4, Padi4, Peptidylarginine deiminase IV, Protein-arginine deiminase type IV

