

## CD40 Recombinant Protein

### OVERVIEW

---

|                           |                      |
|---------------------------|----------------------|
| <b>PRODUCT REFERENCE</b>  | PX-P4082             |
| <b>ORIGIN SPECIES</b>     | Homo sapiens (Human) |
| <b>FRAGMENT</b>           | Met1~Arg193          |
| <b>TAG</b>                | C-terminal His Tag   |
| <b>DELIVERY CONDITION</b> | Dry Ice              |

### PRODUCT INFORMATION

---

|                                |                                                                                                   |
|--------------------------------|---------------------------------------------------------------------------------------------------|
| <b>EXPRESSION SYSTEM</b>       | Eukaryotic expression                                                                             |
| <b>HOST</b>                    | Mammalian cells                                                                                   |
| <b>MOLECULAR WEIGHT</b>        | 22.32kDa                                                                                          |
| <b>FORM</b>                    |                                                                                                   |
| <b>BUFFER</b>                  | PBS pH7.5                                                                                         |
| <b>STABILITY &amp; STORAGE</b> | Store at - 20°C to -80°C.<br>It is recommended that the protein be aliquoted for optimal storage. |
| <b>UNIPROT</b>                 | P25942                                                                                            |
| <b>ALIAS/SYNONYMS</b>          | TNFRSF5                                                                                           |

