



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

CD66c Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P4088
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Gly320
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	35-48kDa
FORM	
BUFFER	PBS pH7.5
STABILITY & STORAGE	Store at - 20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage.
UNIPROT	P40199
ALIAS/SYNONYMS	CD66c,CEACAM6, NCA



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.