



Proteogenix

A n t i b o d i e s a n d P r o t e i n s

Proteogenix SAS
Espace Européen de l'Entreprise
19 rue de La Haye
67300 Schiltigheim, France
Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88
contact@proteogenix.fr
www.Proteogenix.science

CD66c Recombinant Protein

OVERVIEW

CATALOG NO	PX-P4088
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Gly320
MOLECULAR WEIGHT	35-48kDa
DELIVERED WITH TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
PURITY	>90% by SDS-PAGE
UNIPROT ID	P40199
FORM	Liquid
BUFFER	C-terminal His Tag
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)

More info

SYNONYMS

CD66c, CEACAM6, NCA

UNIPROT LINK

<http://www.uniprot.org/uni/prot/P40199>

SEQUENCE

MGPPSAPP CRLHVPWKEVLLTASLLTFWNPPTAKLTIESTPFNVAEGKEVLLLAHNLPQNRIGYSWYKGERVDGNSLIVGYVIGTQQATPG
PAYSGRETIYPNASLLIQNVTDGTYTLQVIKSDLVNEEATGQFHVYPELKPSSNNNSNPVEDKDAVAFTCEPEVQNTTYLWWVNGQS
LPVSPRLQLSNGNMTLTLLSVKRNDAGSYECEIQNPASANRSDPVTLNVLYGPDGPTISPSKANYRPGENLNLSCHAASNPPAQYSWFINGT
FQQSTQELFIPNITVNNSGSYMCQAHNSATGLNRRTVTMITVSG