

Human CX3CR1 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P4093
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1~Pro247
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	47kDa
FORM	
BUFFER	PBS pH7.5
STABILITY & STORAGE	Store at - 20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage.
UNIPROT	P42081
ALIAS/SYNONYMS	CD28LG2

