



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

CD91 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P4094
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Ala20~Gly292
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	37kDa
FORM	
BUFFER	PBS pH7.5
STABILITY & STORAGE	Store at - 20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage.
UNIPROT	Q07954
ALIAS/SYNONYMS	A2MR,APR



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.