



CD137 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P4107
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1~Gln186
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	30kDa
FORM	
BUFFER	PBS pH7.5
STABILITY & STORAGE	Store at - 20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage.
UNIPROT	Q07011
ALIAS/SYNONYMS	CD137, ILA, 4-1BB

