

CD253 Recombinant Protein

OVERVIEW

PRODUCT REFERENCE	PX-P4111
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1~Gly281
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	33.78kDa
FORM	
BUFFER	PBS pH7.5
STABILITY & STORAGE	Store at - 20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage.
UNIPROT	P50591
ALIAS/SYNONYMS	APO2L, TRAIL

