



## CD270 Recombinant Protein

### OVERVIEW

---

|                           |                      |
|---------------------------|----------------------|
| <b>PRODUCT REFERENCE</b>  | PX-P4113             |
| <b>ORIGIN SPECIES</b>     | Homo sapiens (Human) |
| <b>FRAGMENT</b>           | Met1~Val202          |
| <b>TAG</b>                | C-terminal His Tag   |
| <b>DELIVERY CONDITION</b> | Dry Ice              |

### PRODUCT INFORMATION

---

|                                |   |
|--------------------------------|---|
| <b>EXPRESSION SYSTEM</b>       | Eukaryotic expression   |
| <b>HOST</b>                    | Mammalian cells   |
| <b>MOLECULAR WEIGHT</b>        | 36kDa   |
| <b>FORM</b>                    | Frozen  |
| <b>BUFFER</b>                  | PBS, pH7.5  |
| <b>STABILITY &amp; STORAGE</b> | Store at - 20°C to -80°C.<br>It is recommended that the protein be aliquoted for optimal storage. |
| <b>UNIPROT</b>                 | Q92956  |
| <b>ALIAS/SYNONYMS</b>          | HVEA, HVEM,TR2  |

