



## PAPD7 MOUSE Protein

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-P4254
<b>ORIGIN SPECIES</b>	House Mouse
<b>FRAGMENT</b>	Ser2-Arg542
<b>TAG</b>	
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Prokaryotic expression
<b>HOST</b>	Escherichia coli (E.coli)
<b>MOLECULAR WEIGHT</b>	61.05kDa
<b>FORM</b>	
<b>BUFFER</b>	TBS pH7.5 with urea 4M
<b>STABILITY &amp; STORAGE</b>	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
<b>UNIPROT</b>	Q6PB75
<b>ALIAS/SYNONYMS</b>	Terminal nucleotidyltransferase 4A, DNA polymerase sigma, Non-canonical poly(A) RNA polymerase PAPD7, PAP-associated domain-containing protein 7, TRAMP-like complex polyadenylate polymerase, Terminal guanylyltransferase

