

6 kDa early secretory antigenic target(esxA)

OVERVIEW

PRODUCT REFERENCE	PX-P4331
ORIGIN SPECIES	
FRAGMENT	Met1-Gln70
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	33.87kDa
FORM	
BUFFER	TBS pH8.0
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P9WNK7
ALIAS/SYNONYMS	esaT6,ESAT-6

